

EPOND Part-No.	OEM Part-No.	Description
Nebulizers/Spray chambers/accessories		
Quartz Nebulizer		
EP 2050-SILQ-50-05		MEINHARD® SilQ+ Ultra high purity quartz nebulizer
EP 2050-TQP30-A05	N8142098	MEINHARD® TQ+ Ultra high purity quartz nebulizer, 0.5 ml/min uptake, 30 PSI
EP 2050-TQP50-A05	N8122461	MEINHARD® TQ+ Ultra high purity quartz nebulizer, 0.5 ml/min uptake, 50 PSI
EP 2050-MQ30-A05		MEINHARD® Type A concentric quartz nebulizer, 0.5 ml/min uptake, 30 PSI
EP 2050-MQ30-A1		MEINHARD® Type A concentric quartz nebulizer, 1 ml/min uptake
EP 2050-MQ30-A3		MEINHARD® Type A concentric quartz nebulizer, 3 ml/min uptake
EP 2050-MQ50-A05		MEINHARD® Type A concentric quartz nebulizer, 0.5 ml/min uptake, 50 PSI
Glass Nebulizer		
EP 2050-TRP050-05		MEINHARD® TR+ glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-A05		MEINHARD® Type A glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-A1		MEINHARD® Type A glass nebulizer, 1 ml/min uptake
EP 2050-M030-C1	N8102011	MEINHARD® Type C glass nebulizer, 1 ml/min uptake
EP 2050-M030-K1		MEINHARD® Type K glass nebulizer, 1 ml/min uptake
EP 2050-M050-C05	N8145012	MEINHARD® Type C glass nebulizer, 0.5 ml/min uptake
PFA Nebulizers with Integrated Capillaries		
EP 2550-2020		PFA-20 MicroFlow nebulizer, HF resistant, 20µl/min uptake
EP 2550-2002		PFA-100 MicroFlow nebulizer, HF resistant, 100µl/min uptake
EP 2550-2027	N0777376	High temperature APEX-100 PFA nebulizer
EP 2550-MEST384	N8145101	ST3 MicroFlow PFA nebulizer
PFA-ST Nebulizer		
ES-2040		PFA-ST MicroFlow low pressure nebulizer
Poly-Pro Nebulizer		
EP 2550-4040-28		Poly-Pro-ST self-aspirating nebulizer
ES-2040-7000	N0777375	Apex-ST PFA MicroFlow nebulizer
PEEK Nebulizer		
EP 3150-0720		PEEK MiraMist nebulizer
Nebulizer Cleaner		
PRL02		Pearl Nebulizer Cleaner for all Quartz and glass concentric nebulizers
Internal Standard Addition Tee for ESI ST Nebulizers		
EP 0317-1250		0.25mm I.D. internal standard addition tee for ESI ST nebulizers, 14cm
EP 0317-1500		0.5mm I.D. internal standard addition tee for ESI ST nebulizers, 14cm
Nebulizer Gas Quick Connects		
EP-F4		Gas Quick connect for MEINHARD® nebulizers with PFA-Press Fit for 4 mm
EP-F4-54		ICPMS Nebulizer Gas Quick Connect for Nexion/Elan with PFA-Press Fit for 4mm tubing
Cyclonic Spray Chambers		
EP 2040-0200	WE025221	Quartz cyclonic spray chamber for ball joint injector
EP 2040-0200ES	WE025221	Quartz cyclonic spray chamber for ball joint injector with Easy-Seal inlet for PE NexION
EP 2040-0250ESQ	WE8145013	Quartz cyclonic spray chamber for ball joint injector with Easy-Seal inlet for PE NexION
EP 2040-0260ESP	WE8145014	Pyrex cyclonic spray chamber for ball joint injector with Easy-Seal inlet for PE NexION

Torches/shields**Quartz Torches, Demountable Injector**

EP 1102-1000	N8145106	ESI precision o-ring-free quartz torch, demountable, axial
EP 2040-0100	N8122006	Quartz torch, demountable
EP 2040-0120	W1008384	Quartz torch, demountable, high efficiency
Quartz One-Piece Torch with Ball Joint		
EP 2040-0090	N8102001	Quartz torch, fixed injector, long tube

Injectors and accessories**Platinum Injectors**

EP 1113-3100	N8145107	Platinum injectors, 1mm, cassette mount
EP 1113-3200	N8122412	Platinum injectors, 2mm, cassette mount

Quartz Injectors for Standard Torch

EP 2032-0330	N0681631	Quartz injector for standard torch, 1.2mm, bayonet mount
EP 2032-0350	N8125029	Quartz injector for standard torch, 2mm, bayonet mount

Platinum Injectors for O-Ring-Free Torch

EP 1113-2200	N8145111	Platinum injector for o-ring-free torch, 2mm, cassette mount
--------------	----------	--

Alumina Injectors

EP 2040-0186	N8126039	Alumina injector, 0.85mm
EP 2040-0185	N8126040	Alumina injector, 1.5mm
EP 2032-0530	N8126041	Alumina injector, 2mm

Quartz Injectors with Ball Joint

EP 1123-1101	WE027030	Quartz injector with ball joint, 0.85mm, support ring
EP 1123-1151	WE027005	Quartz injector with ball joint, 1.5mm, support ring
EP 1123-1201	WE023948	Quartz injector with ball joint, 2mm, support ring

Quartz Injectors for O-Ring Free Torch

ES-1124-2150	N8145115	Quartz injector for o-ring-free torch, 1.5mm, cassette mount
ES-1124-2250	N8145117	Quartz injector for o-ring-free torch, 2.5mm, cassette mount

Cones**Sampling Cones**

EP 2040-0612	W1033612	Ni sampler Cone, Ni base
EP 2040-0807	WE013087	Ni sampler Cone, Ni base
SC 2011-Ni	WE021140	Ni sampler cone, Ni base
EP 2040-0614	W1033614	Pt sampler cone, Ni base
EP 2040-0250	N8122025	Pt sampler cone, Ni base
SC 2013-Pt	WE027802	Pt sampler cone, Ni base
SC 5011	N8120511	O-ring, brown Viton for sampler(Qty 5)
EP 2040-0148	W1040148	Gaskets, Aluminum for sampler(Qty 3)

Skimmer Cones

EP 2040-0995	W1033995	Al hyper skimmer cone, Al base
EP 2040-0356	W1026356	Ni skimmer cone, Ni base
EP 2040-0000	N8122000	Ni skimmer cone, Ni base
SC 2012-Ni	WE021137	Ni skimmer cone, Ni base
EP 2040-0907	W1026907	Pt skimmer cone, Ni base
EP 2040-0350	N8122035	Pt skimmer cone, Ni base
SC 2014-Pt	WE027803	Pt skimmer cone, Ni base
SC 5012-995	N8120512	O-rings, black Viton for skimmer(Qty 5)
EP 2040-0123	9902123	O-rings, Viton for hyper skimmer(Qty 2)
EP 2040-0737	9919737	Screw, steel for hyper skimmer(Qty 2)
SC6000	WE17142	Insertion/Removal tool

Cone Cleaning Kit & Detergent

Cone Cleaning Kit

Citranox Detergent

EPMNTKIT1

EPCitranox

Accessories

Solution Inlet Adapter

Solution Inlet Adapter for 4 mm standard liquid input nebulizers

Copper Load Coils

Copper load coil for PerkinElmer Elan 6100DRC, DRC II, DRC-e & NexION

Copper load coil for PerkinElmer Elan 500, 5000, 6000, 9001

ELAN Bayonet HF Resistant Sample Intro Kits

Complete HF-resistant sample introduction kit

Complete high purity sample introduction kit

ELAN Cassette HF Resistant Sample Intro Kits

Complete HF-resistant sample introduction kit

Complete high purity sample introduction kit

NexION HF Resistant Sample Intro Kits

Complete HF-resistant sample introduction kit

Complete high purity sample introduction kit

Autosampler Probes

0.8mm I.D. ESI SC-DX *FAST* autosampler probe, carbon fiber support, 150cm

1mm I.D. ESI SC-DX *FAST* autosampler probe, carbon fiber support, 150cm

0.25mm I.D. ESI SC-DX autosampler probe, carbon fiber support, 150cm

0.5mm I.D. ESI SC-DX autosampler probe, carbon fiber support, 150cm

0.25mm I.D. ESI Non-DX autosampler probe, carbon fiber support, 80cm

0.5mm I.D. ESI Non-DX autosampler probe, carbon fiber support, 80cm

0.25mm I.D. PE S10 autosampler probe, carbon fiber support, 100cm

0.5mm I.D. PE S10 autosampler probe, carbon fiber support, 100cm

Electron Multipliers

ETP Electron Multiplier for PE Elan 6000

ETP Electron Multiplier for PE NexION

ETP Electron Multiplier for PE Elan 6X00/9000/DRC

Pipette Tips

100µl PFA PureTip pipette tips, pre-cleaned, pack of 10

1000µl PFA PureTip pipette tips, pre-cleaned, pack of 10

EP-FK3 WE024372

EP 2040-0816 WE021816

EP 2040-0411 WE012411

ES-2170-0472-19 N8122362

ES-2171-0472-19 N8122363

ES-2170-9019 N8122414

ES-2141-9019 N8122415

EP 2170-0472-19 N8142002

EP 2171-0472-19 N8142001

ES-5037-3755-150

ES-5037-3995-150

ES-5037-3255-150

ES-5037-3505-150

ES-5037-3250-080

ES-5037-3500-080

ES-5037-9250-100

ES-5037-9500-100

EPM 14210 N812500

EPM 45000 N8145000

EPM 14217 N8125050

EP 7000-0101 N0777394

EP 7000-1001 N0777395